

Honeywell Specialty Chemicals Seelze GmbH

Postfach 10 02 62, 30918 Seelze

Tel.: +49 (5137) 999-0, Fax: +49 (5137) 999-123

www.honeywell.com

Honeywell Specialty Chemicals Seelze GmbH • Postfach 10 02 62 • 30918 Seelze

Ihre Zeichen

Ihre Nachricht

Unsere Zeichen

Telefon: 05137 999-

Seelze, 13.05.2020

Hydranal CoE

353

HYDRANAL™ NEXTGEN SC Reagents

To Whomsoever It May Concern

This letter confirms that there is no change in the chemical composition of the Hydranal™ titrating agents when sold as HYDRANAL NEXTGEN SC version. The suffix “SC”, representing “Smart Chemical”, added after the original catalog number of the reagents indicates that an RFID tag has been applied to the reagent bottle.

This RFID tag enables certain data concerning the reagent to be read by the Mettler Toledo SmartChemical Reader. The approach in no way alters the composition of the reagents themselves.

Customers wishing to upgrade their Hydranal reagents to Hydranal SC can do so without any corresponding change in the titration method.

For more detailed information please do not hesitate to contact hydranal@honeywell.com.

Kind regards,

Honeywell Specialty Chemicals Seelze GmbH

Thomas Wendt
Supervisor HYDRANAL Center of Excellence



Dr. Burkhard Vehre
Quality Manager